

Sample preparation guide for sorting

- Test your antibody and fluorophore combination before booking a sort.
- Cells must be in a single cell suspension.
- Filter the sample just prior to sorting in order to remove cell clumps. Use the 5 ml Falcon polystyrene round-bottom tubes with cell-strainer cap (cat# 352235, 35µm nylon mesh) for this purpose.
- Bring the cells in Mg/Ca⁺⁺ free buffer, containing 0.1-1% of BSA or 1-2% of dialyzed FBS. In order to prevent cell clumping during the sort, add EDTA (up to 5mM, prevents cation dependent cell to cell adhesion) and/or DNase (up to 50µg/ml, prevents cell clumping due to a cell death). The cells can also be brought in Phenol Red free media.
- Use Dead Cell Exclusion dyes to eliminate dead cell.
- Sample cell concentration is cell type and nozzle size dependent - recommended concentration 5 -10x10⁶cells/ml.
- Cells can be collected into:
 - 5ml (12x75mm) and 15 ml tubes
 - 6, 24 and 96 well plate
 - Slides
- Collection tube size/ number: a rule of thumb is that with 70µm nozzle we can put up to 3 million droplets in the 12x75mm tube. With 100µm nozzle we can put up to 1 million droplets in 12x75mm tube.
- In order to get a better recovery, pre-coat the collection tubes by filling them with media+serum or serum only for at least 30' prior to sorting.
- Solution to sort into:
 - PBS based buffer containing FBS or BSA
 - Complete media for cell growth with 10-20% FBS
 - TRIzol reagent (for RNA extraction)
 - Etc.
- Provide compensation controls and gating controls.